



Summer

Summer[®]



Table of Contents

Weighting	2
-----------------	---

Weighting



Balança Sartorius TE212

Model: 9811

Built-in application programs: Net-total formulation, weighing in percent, counting, mass unit conversion by toggling, dynamic weighing.

Electronic quadrant scale for samples determining the weight of fabrics for samples of 100 square centimeters, weighing Range 1 -12,000 g/m² (121 g effective) scale 1 g/m² (0.01 g effective) bowl 116 mm diameter.

Balance BS945-15

Model: 9806

Control and industrial balance type BS 945, model BS945-15, 0-15kg / 1 gr. Dimensions of the plate 215 x216 mm.



Cutting blades

Model: 9815/2

Cutting blades for sample cutter 9815 and 9815/1. Box with 100 pieces.



Sample cutter

Model: 9815

Sample cutter for precise cutting of circular samples of exactly 100 square centimeters, for use on quadrant scales. Accurate cutting edges by circular guided blades, and by draw-cut. 4 adjustable and interchangeable blades, with pad of porous rubber.



